

ONCOTYPE MAP™ PAN-CANCER TISSUE TEST

RAPID COMPREHENSIVE GENOMIC PROFILING TO AID THERAPY SELECTION

Every day matters. Why wait?

Oncotype MAP Pan-Cancer Tissue test delivers rapid, comprehensive tumor profiling to aid therapy selection for patients with advanced, metastatic, refractory, relapsed or recurrent cancer. The test identifies actionable genomic alterations within 3-5 business days** to guide timely treatment decisions, usually before the next follow-up visit.

Utilizing next generation sequencing (NGS) and a broad array of immunohistochemical (IHC) stains and panels, the Oncotype MAP Pan-Cancer Tissue test identifies clinically-relevant genomic alterations from minimal amounts of tissue. The NGS test can be run on patient specimens as small as 3 mm² of tissue with ≥15% tumor content (20% optimal for MSI), allowing you to identify appropriate patient therapy based on tumor characteristics.

Results are provided in an easy-to-interpret report with NCCN® actionable indications for therapy, along with potential evidence-based therapies and clinical trials, to guide treatment options for a breadth of solid tumor types.



Answers you need at the speed your patients deserve



Fast

- Results reported in 3-5 business days
- Quantity Not Sufficient (QNS) results typically reported within 48 hours



Comprehensive

- 257 genes, 11 select fusions
- 21 IHC stains including PD-L1 (22C3, SP142)
- 30+ tumor-specific panels
- Copy number variants, tumor mutation burden, microsatellite instability, and protein expression



Actionable

- 100+ therapies, 45+ combination therapies, 650+ clinical trials
- NCCN Compendium recommendations



Accurate

- Leverages a proprietary method designed for small specimens (3 mm² with tissue with ≥15% tumor content), enabling a low failure rate on both FNA and FFPE samples*
- >99% sensitivity for SNVs and Indels, detection down to ≥7.5% mutant allele frequency (MAF) with >99% specificity*

Oncotype MAP Pan-Cancer Tissue test helps guide therapy selection for today and captures emerging markers that may find their way into guidelines tomorrow.

**Order now for your patients
with solid tumors**

Oncotype MAP™ Pan-Cancer Tissue Test Markers

Immunohistochemistry												
Single IHC Stains		Tumor-Specific Panels										
ALK	PD-L1 (22C3)	Anal Carcinoma: PD-L1 (22C3), PD-1, TS, TUBB3 Appendix: HER2, PD-L1 (22C3), TOPI, PTEN Bladder: PD-L1 (22C3), PD-L1 (SP142), hENT1, TUBB3 Bone Cancer: TOPI, MGMT, CAIX, HENT1 Breast: AR, PD-L1(SPI42), PD-L1 (22C3), TP CNS/Brain Cancers: MGMT, CAIX, TUBB3, TOPI Cervical: PD-L1 (22C3), CAIX, HENT1, TOPI Colorectal: HER2, PD-L1 (22C3), TOPI, PTEN, TS Gastric/Esophageal: HER2, PD-L1 (22C3), PTEN, TS GIST: PD-L1 (22C3), MET Head and Neck: PD-L1 (22C3), CAIX, TUBB3, PTEN Head and Neck Salivary Gland: HER2, AR, CAIX, PTEN Hepatobiliary/Cholangiocarcinoma: hENT1, HER2, TP, PD-L1 (22C3) Hepatobiliary/Gallbladder: hENT1, HER2, PD-L1 (22C3), TOPI Hepatobiliary/Hepatocellular: hENT1, PD-L1 (22C3), CAIX, MET Kidney: PD-L1 (22C3), MET, CAIX, HENT1 Melanoma: PD-L1 (22C3), MGMT, PTEN, TUBB3 Mesothelioma: PD-L1 (22C3), TS, hENT1, TUBB3						MMR: MLH1, MSH2, MSH6, PMS2 Neuroendocrine: PD-L1 (22C3), MGMT, PTEN, TP NSCLC: PD-L1 (22C3), ALK, PD-L1 (SP142), PTEN, TS Ovarian: ER, HER2, TOPI, TUBB3 Pancreatic: hENT1, PTEN, TP, TOPI Penile Cancer: PD-L1 (22C3), TP, CAIX, TUBB3 Prostate: AR, PTEN, TUBB3, PD-L1 (22C3) Sarcoma: CAIX, TUBB3, TOPI, MGMT SCLC: PD-L1 (22C3), TOPI, MGMT Small Bowel: TOPI, CAIX, TUBB3, hENT1 Testicular Cancer: PD-L1 (22C3), TUBB3, hENT1 Thyroid: PD-L1 (22C3), ALK, TUBB3, CAIX Thymoma/Thymic Carcinoma: PD-L1 (22C3), TUBB3, hENT1, TS Uterine: PD-L1 (22C3), HER2, MGMT, TUBB3 Uveal Melanoma: PD-L1 (22C3), MGMT, TUBB3 Vulvar Cancer: PD-L1 (22C3), CAIX, TUBB3, hENT1 Other Solid Tumors: PD-L1 (22C3), HER2, TOPI, PTEN, TS				
AR	PD-L1 (SP142)											
CAIX	PR											
ER	PTEN											
hENT1	ROS1											
HER2	TOPI											
IDO	TP											
MET	TRKpan											
MGMT	TS											
PD1	TUBB3											
MMR Panel												
MLH1	MSH6											
MSH2	PMS2											
257 Gene NGS Panel												
ABC1	AREG	BUB1B	CYP19A1	EWSR1	GAS6	KEAP1	MYC	PIK3R1	RPTOR	TGFB2		
ABCC1	ARID1A	CALR	CYP11A1	EZH2	GATA3	KIT	MYCN	PIM1	RRM1	TGFB3		
ABCC2	ARID1B	CBL	CYP2D6	FAM175A	GLI1	KRAS	MYOD1	PLCB4	SDHB	TGFBR1		
ABL1	ARID2	CCND1	CYP3A4	FANCA	GNA11	MAF	NBN	PLCG1	SDHC	TGFBR2		
ACVR1	ATM	CCND2	CYSLTR2	FANCC	GNAQ	MAP2K1	NF1	PMS2	SETD2	TNFAIP3		
ACVR1B	ATR	CCND3	DCK	FANCD2	GNAS	MAP2K2	NF2	POLD1	SF3B1	TNK1		
ACVR2A	ATRX	CCNE1	DDR2	FANCE	GSTP1	MAP3K1	NFE2L2	POLE	SMAD1	TOP2A		
ACVR2B	AURKA	CD274	DICER1	FANCF	HAMP	MAPK1	NOTCH1	PP2R1A	SMAD2	TP53		
ACVRL1	AURKB	CDA	DNMT3A	FANCG	HDAC2	MAPK3	NOTCH2	PTCH1	SMAD4	TSC1		
ADAMTS1	AXIN1	CDC73	EGFR	FANCM	HGF	MAPKAPK5	NOTCH3	PTEN	SMAD5	TSC2		
ADAMTS6	AXL	CDH1	EMSY	FAT1	HNF1A	MDM2	NPM1	PTPN11	SMAD9	TSHR		
ADAMTS9	B2M	CDK4	EP300	FBXW7	HRAS	MDM4	NRAS	RAD50	SMARCA4	TYMS		
ADAMTS16	BAP1	CKD6	EPCAM	FCGR2A	HSD3B1	MED12	NTRK1	RAD51C	SMARCB1	VEGFA		
ADAMTS18	BAR1	CDK12	EPHA5	FGD4	IDH1	MEN1	NTRK2	RAD51D	SMO	VHL		
ADAMTSL1	BCOR	CDKN2A	EPHA7	FGF3	IDH2	MET	NTRK3	RAF1	SOCS1	WT1		
AKT1	BMP6	CHEK1	ERBB2	FGF4	IGF1R	MGMT	PALB2	RB1	SPOP	XRCC1		
AKT2	BMPR1A	CHEK2	ERBB3	FGFR1	IKZF1	MLH1	PBRM1	RBM10	STAG2	YES1		
AKT3	BMPR1B	CHFR	ERBB4	FGFR2	IL6R	MPL	PDCD1LG2	RECQL	STAT3			
ALK	BNIP3	CHKA	ERCC1	FGFR3	JAK1	MRE11A	PDGFRA	RET	STAT5A			
AMER1	BRAF	CIC	ERCC2	FGFR4	JAK2	MSH2	PDGFRB	RHEB	STAT5B			
APC	BRCA1	CREBBP	ERCC3	FLT3	JAK3	MSH6	PIK3CA	RICTOR	STK11			
APLNLR	BRCA2	CSF1R	ERRF1	FLT4	KDM5C	MTHFR	PIK3CB	RIT1	SUFU			
AR	BRIP1	CTLA4	ESR1	FOXL2	KDM6A	MTOR	PIK3CD	RNF43	TERT-p			
ARAF	BTK	CTNBN1	ESR2	FUBP1	KDR	MUTYH	PIK3CG	ROS1	TGFB1			

Genetic Structures Tested: Single nucleotide variants (SNVs) and insertions/deletions in coding regions of genes listed above; UTRs and splice junctions when actionable (e.g. MET exon 14 skipping and EGFRvIII). MSI; mutation burden (SNV's, insertions, deletions) based on -1 megabase; select fusions involving ALK, BRAF, FGFR1, FGFR2, FGFR3, MET, RET, ROS1, NTRK1, NTRK2, NTRK3 (ETV6); and copy number variants.

Turnaround time is based on when qualified sample is received. Mutation calls may not be available from some regions due to pseudogenes or sequence context. Select IHCs may not be run if already performed within the last six months unless indicated in the notes section. HER2 equivocal by IHC will be reflexed to FISH testing in select tumor types. Reflex testing will exceed standard turnaround time for results. MMR includes the following IHCs: MLH1, MSH2, MSH6, PMS2.



Clinical Laboratory Improvement Amendments



Contact Paradigm Diagnostics/Exact Sciences
 1-866-662-6897 | oncotypemap@exactsciences.com
www.oncotypemap.com

About Exact Sciences

A leading provider of cancer screening and diagnostic tests, Exact Sciences helps people get the answers they need to make more informed decisions across the cancer continuum. Building on the success of the Cologuard® and Oncotype DX® tests, Exact Sciences is investing in its product pipeline to take on some of the deadliest cancers and improve patient care. Through an innovative, rigorous approach, and with the support of visionary collaborators, we're helping advance the fight against cancer.

References:

* Data on file, in-house assay optimization protected by trade-secret/patent regulations.
 † Turnaround time is based on qualified sample receipt.
 ‡ Morris S, Subramanian J, Gel E, Rungger G, Thompson E, Mallery D, et al. Performance of next-generation sequencing on small tumor specimens and/or low tumor content samples using a commercially available platform. PLoS ONE. (2018); 13(4): e0196556. <https://doi.org/10.1371/journal.pone.0196556>.

NCCN and NCCN Guidelines are trademarks of the National Comprehensive Cancer Network. NCCN makes no warranties of any kind whatsoever regarding their content, use or application and disclaims any responsibility for their application or use in any way.

Oncotype MAP and Oncotype MAP Pan-Cancer Tissue Test are trademarks of Genomic Health, Inc. Exact Sciences is a registered trademark of Exact Sciences Corporation. ©2020 Genomic Health, Inc. All rights reserved. EXS50041_0121